

Website Layout & Copy

Get your beautiful, healthy smiles...faster!

AcceleDent® Optima™ gives patients faster and more comfortable orthodontic treatment with braces or aligners. AcceleDent Optima is the most advanced accelerated orthodontic device with a comprehensive suite of technology driven features. Simple and easy-to-use, AcceleDent Optima is clinically proven to speed up tooth movement while reducing orthodontic discomfort.



AcceleDent®
OPTIMA™

Healthy, Beautiful Smiles...Faster™!

What Is AcceleDent Optima?

AcceleDent Optima is clinically proven to accelerate tooth movement by up to 50%. It's the first orthodontic device of its kind to offer accelerated orthodontic treatment and help decrease orthodontic discomfort in a noninvasive, nonsurgical way. AcceleDent Optima is also the first and only accelerated orthodontic device to connect patients and practices with an app to monitor, monitoring, direct messaging & rewards. AcceleDent is only available by prescription from your doctor.

Safe & Effective

Leading doctors around the world use and endorse AcceleDent Optima's SoftPulse Technology® as safe and effective in accelerating tooth movement and reducing the discomfort associated with orthodontics.

Lightweight & Easy-to-use

AcceleDent Optima was designed to fit your lifestyle. Use the hands-free device for just 20 minutes a day during routine activities such as reading, watching TV, listening to music, browsing the internet or showering.



[Link to Introducing AcceleDent Optima Video - <https://youtu.be/Hn9xLQYvdi8>]

To embed the video on your website, copy and paste the HTML code below:

```
<iframe width="560" height="315" src="https://www.youtube.com/embed/Hn9xLQYvdi8" frameborder="0" allow="autoplay; encrypted-media" allowfullscreen></iframe>
```

How Does AcceleDent Optima Work?

By inserting the mouthpiece fitted around your braces or aligners for just 20 minutes a day, AcceleDent Optima speeds tooth movement through the use of SoftPulse Technology®. AcceleDent Optima's patented SoftPulse Technology applies precisely calibrated vibrations called micropulses that transmit through the roots of your teeth to the surrounding bone. This gentle vibration helps increase the cellular activity and speeds the rate in which your teeth can move. AcceleDent's gentle micropulses can also reduce orthodontic discomfort for patients of all ages, including children, teens and adults.

Designed to Fit Your Lifestyle

AcceleDent Optima combines high-tech performance features and a comprehensive app to enhance your accelerated orthodontic experience and empower you to become a full participant in your treatment.



Features

- Small & Lightweight Design
- Waterproof Activator
- Wireless Charging
- Bluetooth® Enabled

AcceleDent App

- Dashboard reports for tracking progress (also accessible on MyAcceleDent.com)
- Motivation tools to encourage usage
- Community connectivity to share progress and updates with friends, family and the treatment team

What are AcceleDent Patients Saying?

- “AcceleDent was amazing for me because I wanted to accelerate the treatment process and it was ideal to reduce the discomfort I was experiencing.” – Karen M.
- “Thanks to AcceleDent, what should have been a two-year ordeal took only one year. I couldn’t recommend AcceleDent more!” – Miguel A.
- “I have used AcceleDent every day during my orthodontic treatment, and I love it!” – Jennifer K.
- “I am so thankful to have found AcceleDent. I especially loved using it on days when I would switch trays, which helped alleviate the pain. I am so impressed with how quickly I saw results, and would absolutely recommend it to everyone!” – Andrea A.
- “AcceleDent is my best friend on adjustment day!” – Sonyah M., Patient
- “Not only have my teeth moved very quickly, but using AcceleDent has greatly reduced the discomfort that I feel while wearing braces.” – Barbara K.
- “I’m getting amazing results with AcceleDent!” – Scott M.
- “I strongly recommend AcceleDent for anyone who wants to reduce their treatment time, and it’s so easy to use!” – Jean B.

How to Use and Care for AcceleDent Optima and Use the AcceleDent App



[Link to AcceleDent Optima Patient Video - https://youtu.be/vE_ujubZrF4]

To embed the video on your website, copy and paste the HTML code below:

```
<iframe width="560" height="315" src="https://www.youtube.com/embed/vE_ujubZrF4" frameborder="0" allow="autoplay; encrypted-media" allowfullscreen></iframe>
```